

DETERMINATION OF ORGANOCHLORINE PESTICIDE IN EDIBLE OIL (USING SWEEP CO-DISTILLATION CLEAN-UP METHOD)

by: HALIMAH, M; M D PAUZI, A; SOH, S C and AINIE, K



295

MPOB INFORMATION SERIES • ISSN 1511-7871 • JUNE 2005

MPOB TT No. 298

This method is used to determine the concentration of 15 organochlorine pesticides in edible oil using capillary column with electron detector (ECD) (Figure 1).



Figure 1. Unitrex for organochlorine analysis.

DEFINITION

The compounds as shown in Table 1 are to be determined using this method.

PRINCIPLE

The method involves spreading a thin film of the sample over glass surfaces (glass beads)

TABLE 1.

Compound	CAS registry No.
α -BHC	319-84-6
β -BHC	319-85-7
γ -BHC (lindane)	58-89-9
δ -BHC	319-86-8
α -Chlordane	5103-71-9
γ -Chlordane	5103-74-2
α , p , p' -DDD	72-54-8
p , p' -DDE	72-55-9
p , p' -DDT	50-29-3
Endrin	72-20-8
Endrin ketone	53494-70-5
Heptachlor	76-44-8
Heptachlor epoxide	1024-57-3
Methoxychlor	72-43-5

inside a fractionation tube, the pesticides are then partitioned into the gas stream and are carried up and out of the sidearm of the fractionation tube into a trap packed with sodium sulphate and partially deactivated Florisil. The pesticides are then eluted from the trap and determined by GC-ECD (Figure 2).




Figure 2. Gas chromatogram of 15 organochlorine pesticide reference standards: 1 μ l injection volume of 1 μ g ml⁻¹ concentration. Peaks: IS1 = tetrachloro-*m*-xylene; 1= α -BHC = β -BHC; 3 = lindane; 4 = δ BHC; 5 = heptachlor; 6 = heptachlor epoxide; 7 = γ -chlordane; 9 = p , p' -DDE; 10 = dieldrin; 11 = endrin; 12 = p , p' -DDD, 13 = p , p' -DDT; 14 = endrin ketone; 15 = methoxychlor; IS2 = decachlorobiphenyl.



REFERENCES

HALIMAH, M; M D PAUZI, A; TAN, Y A and CHIAN, S S (2004). Optimization of the sweep co-distillation clean-up method for the determination of organochlorine pesticide residues in palm oil. *J. Oil Palm Research Vol. 16 No. 2*: 30-36.

LUKE, B G; RICHARDS, J C and DAWES, E F (1984). Recent advances in clean-up of fats by

sweep co-distillation. Part 1: organochlorine residues. *J. Assoc. Official Anal. Chem.*, 67: 30-36.

M D PAUZI, A ; HALIMAH, M; AINIE, K and CHIAN, S S (2003). Comparison of solid phase extraction and sweep co-distillation techniques for the determination of organochlorine pesticide residues in refined palm oil. *Sains Malaysiana*, 32(2): 131-142.

For more information kindly contact:

Director-General
MPOB

P. O. Box 10620

50720 Kuala Lumpur, Malaysia.

Tel: 03-89259155, 89259775

Website: <http://mpob.gov.my>

Telefax: 03-89259446