

# DETERMINATION OF HYDROCARBONS (n-alkanes) IN OIL MATRIX

by: NORIZAH HALIM and AINIE KUNTOM



MPOB INFORMATION SERIES • ISSN 1511-7871 • JUNE 2005

MPOB TT No. 300

## SCOPE

This test method describes a method for the determination of hydrocarbons (n-alkanes) in vegetable oils and fats using gas chromatography-mass spectrometric (GCMS) detection. Vegetable oils may be contaminated by mineral fuel oils (diesel-oil) due to the use of transport/storage facilities which are not strictly for transport/storage of food-grade, or due to other reason. The mineral fuel oil consists mainly of hydrocarbons (n-alkanes) with a chain length in between C10 to C24 and their quantification limit is 1.0 mg kg<sup>-1</sup>.

## DEFINITION

Alkanes are a homologous series of saturated aliphatic hydrocarbons. Alkanes have the general formula C<sub>n</sub>H<sub>2n+2</sub> (where n is an integer greater than or equal to 1) and are generally grouped as low molecular weight, e.g. methane, ethane and propane are gases; intermediate molecular weight, e.g. hexane, heptane and acetone are liquids; and high molecular weight, e.g., eicosane (C<sub>20</sub>H<sub>42</sub>) and polyethylene are solids. Paraffin is a mixture of high molecular weight alkanes sometimes called the paraffin series.

## PRINCIPLE

The test sample dissolved in n-pentane is added onto an aluminium oxide (alumina) column, from which the hydrocarbon are eluted with n-pentane (Figure 1). The eluate is analysed using GC-MS (gas chromatography mass spectrometry) (Figures 2 and 3). Comparison is made with a sample of a standard diesel oil.




Figure 1. Set-up for separation of the hydrocarbons from oil matrix.



Figure 2. GC-MSD for analysis of the hydrocarbons.

ISSN 1511-7871



9 771511 787001






Figure 3. Chromatogram of standard n-alkanes.

## RECOVERY

Recoveries of the hydrocarbons from spiked oil matrix at the range of 25-100  $\mu\text{g ml}^{-1}$  were 87.1%-103.5%.

Coefficient of variation was <10% at all concentrations.

Limit of detection was 1  $\mu\text{g ml}^{-1}$ .

For more information kindly contact:

Director-General  
MPOB

P. O. Box 10620

50720 Kuala Lumpur, Malaysia.

Tel: 03-89259155, 89259775

Website: <http://mpob.gov.my>

Telefax: 03-89259446