

# OIL PALM LIQUID CULTURE - MPOB PROTOCOL

by: **AHMAD TARMIZI HASHIM**

MAY 2002

**155**

**MPOB TT No. 138**

MPOB INFORMATION SERIES

ISSN 1511-7871


**A**s an alternative to culture on solid medium, the use of liquid culture for oil palm propagation has been attempted in order to explore its potential for automation and scaling-up propagule production as well as for improving growth and culture uniformity (de Touchet *et al.*, 1991; Teixeira *et al.*, 1995; Wong *et al.*, 1999). However, a detailed protocol of the method has not been published to date. Therefore, the purpose of this technology transfer is the provision of the established MPOB protocol to interested parties.

## OIL PALM LIQUID CULTURE PROTOCOL

### Materials

#### 1. Selection of suitable callus.

Friable embryogenic callus is the best starting inoculum (*Figure 1*).



*Figure 1. Friable embryogenic callus.*

#### 2. Media preparation.

##### Medium 1:

Liquid medium - for culture establishment and maintenance.

Basal medium: MS (Murashige and Skoog, 1962) macro, MS micro, Vitamins, 0.1 g litre<sup>-1</sup> myo-inositol, 0.1 g litre<sup>-1</sup> L-glutamine, 3% sucrose (castor sugar brand MSM).

Hormone: 1 mg litre<sup>-1</sup> 2,4-D or 1 mg litre<sup>-1</sup> 2,4-D + 0.1 mg litre<sup>-1</sup> NAA.

Volume: 20 ml in 100 ml flask.

##### Medium 2:

Liquid induction medium - for maturation induction of culture aggregates.

Basal medium: MS macro, MS micro, vitamins, 0.1 g litre<sup>-1</sup> myo-inositol, 0.1 g litre<sup>-1</sup> L-glutamine, 3% sucrose (castor sugar brand MSM).

No hormone.

#### 3. Preparation of apparatus for sterile work.

Apparatus required are:

- polypropylene sheets (10 x 10.5 cm);
- beakers;
- strainers (sieves) 2 mm, 1 mm, 500 µm, 106 µm mesh sizes;
- forceps/scalpels/spatulas/scoops; and
- all apparatus need to be autoclaved.



## Methodology

### 1. Inoculation.

- a) Transfer 0.03-0.09 g of friable embryogenic callus into 20 ml of MS media.
- b) Replace cotton plugs with polypropylene sheets
- c) Transfer flasks to shakers (*Figure 2*).


Physical conditions: dark room, temperature 27°C ± 1°C. Speed: 90 -100 rpm.



*Figure 2. Flasks on shakers.*

### 2. Subculturing at monthly intervals.

Duration for establishment of liquid cultures may vary between clones (usually more than three subcultures) (*Figure 3*).



*Figure 3. Established cell aggregates.*

- Select cultures (various sizes);
- Streak to check contamination;
- Sieve according to sizes, e.g. 2 mm, 1 mm, 500 µm; and
- Transfer 0.3-1.0 g of cell aggregates to 20 ml MS medium.

### 3. Maturation.

- Transfer to basal liquid MS medium to induce maturation before regeneration (one month).

### 4. Regeneration.

- Transfer cultures to basal MS solid medium (petri dishes, baby jars or flasks).
5. Multiplication of embryoids/shooting/rooting and establishment in soil will follow standard protocol as reported for solid culture system (Zamzuri *et al.*, 1999).

## CONCLUSION

This protocol has been established and used routinely at MPOB. However, the performance of the cultures is genotype dependent and the duration for establishment of liquid cultures varies from clone to clone.

## REFERENCES

- de TOUCHET, B; DUVAL, Y and PANNETIER, C (1991). Plant regeneration from embryogenic suspension culture of oil palm (*Elaeis guineensis* Jacq.). *Plant Cell Report*, 10:529-532.
- MURASHIGE, T and SKOOG, F (1962). A revised medium for rapid growth and bioassays with tobacco tissue cultures. *Physiologia Plantarum*, 15: 473-497.
- TEIXEIRA, J B; SONDAHL, M R; NAKAMURA, T and KIRBY, E G (1995). Establishment of oil palm suspensions and plant regeneration. *Plant Cell, Tissue and Organ Culture*, 40: 105-111.
- WONG, G; CHONG, S P; TAN, C C and SOH, A C (1999). Liquid suspension culture - a potential technique for mass production of oil palm clones. *Proc. of the 1999 PORIM International Palm Oil Congress - Emerging Technologies and Opportunities in the Next Millennium*. Kuala Lumpur. p. 3-11.
- ZAMZURI, I; MOHD ARIF, S; RAJANAIDU, N and ROHANI, O (1999). Commercial feasibility of clonal oil palm planting material production. *PORIM Occasional Paper No. 40*: 52 pp.

For more information kindly contact:

Director-General  
MPOB  
P. O. Box 10620  
50720 Kuala Lumpur, Malaysia.  
Tel: 03-89259155, 89259775,  
Homepage: <http://mpob.gov.my>  
Telefax: 03-89259446