

LASER SPECTROFLUORIMETER SYSTEM

by: TAN YEW AI, CHONG CHIEW LET
LEE CHAK KHIAM¹, HARLINA DAMIRI¹ AND LOW KUM SANG²

JUNE 1999

76

PORIM TT No. 51


PORIM INFORMATION SERIES

ISSN 0128-5726

The *Laser Induced Fluorescence Spectrofluorimeter (LIFmeter) is designed to detect chlorophyll content in both crude and refined vegetable oils. It exhibits a very high sensitivity compared to the

conventional spectrophotometer method.

Operation is as simple as filling a cuvette with the test sample, inserting it into the cuvette chamber and pressing the measure button.



* Application for Utility Innovation pending (Malaysia)

ISSN 0128-5726



Palm Oil Research Institute of Malaysia, Ministry of Primary Industries, Malaysia
P. O. Box 10620, 50720 Kuala Lumpur, Malaysia. Tel: 03-8259155, 8259775, Homepage: <http://porim.gov.my> Telefax: 03-8259446

- 1 Fotontec Sdn. Bhd.
- 2 Universiti Malaya



FEATURES

- Simple drop in measurement
- Compact and portable
- Minimum detection level of 10µg/kg
- High throughput, less then 5 seconds per measurement
- Menu guided operation
- Ultra high linearity throughout measurement range achieved through the use of proprietary modified dual slope measurement technique

SPECIFICATIONS

Light Source	5mW diode laser
Detector	Photo Diode
Detection Limit	10µg/kg of chlorophyll
Display	LCD
Measurement Time	<5 seconds
Dimensions (WxLxH)	220mm x 170mm x 80mm
Power Supply	120/230 V
Weight	1kg

Specifications and design are subject to change without notice

For more information kindly contact:

Director-General
PORIM

P.O. Box 10620
50720 Kuala Lumpur, Malaysia.

or

Fotontec Sdn. Bhd.

No. 2, UM-MTDC Technology Innovation Centre,
University of Malaya, Lembah Pantai
50603 Kuala Lumpur, Malaysia

Telephone: 603-7541078

Facsimile: 603-7541079

E-mail: leec@pc.jaring.my